

**FirstView<sup>®</sup>**


# Rapid test for HIV 1.2 & Syphilis Combo

**For Qualitative  
Detection of HIV 1.2 &  
Syphilis Combo  
Human Serum, Plasma  
or Whole Blood  
Rapid Test**


- In-Vitro diagnostic immune-chromatographic assay for detection of antibodies specific to HIV1/2 and syphilis antibodies (IgM/IgG) to Treponema pallidum.
- Run with Human Serum, Plasma or Whole Blood
- Semi-quantitative point-of-care Test for HIV 1.2 & Syphilis Combo Rapid Test.






## READING TEST RESULTS



Add 1 drops (25 µL) of serum or plasma in sample well "S" using the dropper provided.



Add 2 drop of the Assay buffer



Do not read the results after 30 minutes

ISO 9001  
CERTIFIED

GMP  
CERTIFIED

ISO 13485  
CERTIFIED



**PCTSI**

Imported and Distributed By :

**Phil Capital Trading Solutions Inc.**

904 Atlanta Centre, 31 Annapolis St. Greenhills,  
San Juan City, Metro Manila  
T. (02)7001 0770 | E. philtradinginc6@gmail.com