

FirstView®

One Step Rapid Test Device



Negative
No Infection



IgM Positive
Early primary infection



IgG/IgM Positive
Active primary, repeat infection



IgG Positive
Late stage or, latent infection



Invalid
Re-Test

**Qualitative
Proprietary
Whole Blood
Rapid Test**



Imported and Distributed By :

Phil Capital Trading Solutions Inc.

904 Atlanta Centre, 31 Annapolis St. Greenhills,
San Juan City, Metro Manila

T. (02)7001 0770 | E. philtradinginc6@gmail.com