

U-220

Urine Analyser



Reliable Testing, Intelligent Design

- User-friendly LCD display and menu-driven interface
- Bright, cold light source for clear results
- Long-lasting performance
- Highly sensitive light sensor with interference resistance- Intelligent pH ratio correction
- Fast, quiet thermal printing
- Seamless computer connectivity via RS-232- Storage for 10,000 test results
- Customizable alert thresholds- Adjustable test sensitivity.
- Easily connect external printers for more printing options.
- Strip Available 10,11 Parameter
- Test speed is 120 test 60 test / hour and manually triggered.

ISO 9001
CERTIFIED

GMP
CERTIFIED

ISO 13485
CERTIFIED

U-220 Urine Analyser



SPECIFICATIONS

Test Items:	Urobilinogen, Bilirubin, ketone, Blood, Protein, Micro Albumin, Nitrite, Leucocytes, Glucose, Specific Gravity, pH,
Test Color Wave	470nm, 540nm, 620nm
Test Principle	Reflectance and color analysis
Strips Available	10, 11, Parameters.
Test Speed	120 Tests, 60 Tests/hour and Manually triggered
Display	192 x 64 Graphic LCD display
Data Base Memo	10000 test reports
Printing	Internal 57nm thermal printer
Interface	Standard RS232 interface
USB	Can be connected to the scanner;
Working Environment	Temperature 15-30°C, humidity < 70%
Power Supply	AC220V/110V±15%, 50/60HZ
Power consumption	35VA Size: 253x237x171
Weight	3kg



CE

PCTSI

Imported and Distributed By:

Phil Capital Trading Solutions Inc.

B11 L2 Dove St., Ecoland Subdivision Phase 2, Brgy. Bucana, Davao City, Davao Del Sur

Manila office: 904 Atlanta Centre, 31 Annapolis St.

Greenhills, San Juan City, Metro Manila

T. (02)7001 0770 | E. philtradinginc6@gmail.com