

FirstView®

One Step Rapid Test Device

Bioline
DIAGNOSTICS
Your safety. Our commitment.

Rapid test for Dengue NS1 Ag

**For Qualitative
Detection of
Dengue NS1 Ag
In human serum,
plasma or
whole blood
Rapid Test**

- Dengue NS1 Antigen Test provides an excellent methodology for specifically detecting Dengue NS1 antigen within up to 1 day of infection.
- Run with WHOLE BLOOD or PLASMA or SERUM
- Rapid test result in just 15 -20 minutes
- Semi-quantitative point-of-care Test .




ISO 9001
CERTIFIED



ISO 13485
CERTIFIED

READING TEST RESULTS



Add 3 drops (75µL) of Serum, Plasma, Whole Blood into sample well "S" using the dropper provided.



Add 1 drop (40µL) of the Assay buffer




Result



Positive



Negative



Invalid

PCTSI

Imported and Distributed By:

Phil Capital Trading Solutions Inc.

904 Atlanta Centre, 31 Annapolis St. Greenhills,
San Juan City, Metro Manila

T. (02)7001 0770 | E. philtradinginc6@gmail.com

IVD